
PATIENT AND CAREGIVER USER MANUAL

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www.molnlycke.com/symbols

STERILE EO The device is sterilized using ethylene oxide

Do not reuse

Caution, see instructions for use

Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 13.

Follow instructions for use

Temperature limitation

Keep dry
Keep away from rain

Keep away from sunlight
Keep away from heat

Catalogue number

Use by date / Expiry date

Batch code

Serial number

Single sterile barrier system

Medical Device

Do not use if package is damaged

Manufacturer

MRI unsafe

Humidity limitation

Atmospheric pressure limitation

Applied part type BF

Ingress Protection

ETL Listed Marking

Separate collection of Waste Electrical and Electronic Equipment (WEEE)

Leakage

Blockage


Low battery

4. Precautions

- The Avance Solo Adapt Pump has both visual and audible notifications and alarms. Carry or position the pump so that you will be able to detect audible and visual notifications or alarms.
- Regularly monitor the canister mounted on the pump. If the canister appears full or the pump alarms for blockage, change the canister according to instructions provided in this manual.
- When the Avance Solo Adapt Pump alarms due to a low battery, replace the batteries in the pump according to the instructions provided in this manual. Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care.
- Ensure that the battery lid on the Avance Solo Adapt Pump is closed during therapy.
- If you see signs of possible infection such as pain, reddening, odor, or sensitization of the wound area or sudden change of wound fluid volume or color, immediately contact a healthcare professional.
- Inform your healthcare professional if you are hypersensitive to the product's materials.
- The Avance Solo Canister is provided sterile. Do not use the canister if the inner packaging of the canister is damaged.
- Do not place the Avance Solo Adapt Pump with canister in water or other liquids. If the pump is wet, disconnect the pump and canister and contact your healthcare professional.
- Do not expose the dressing to extensive contact with water. If the transfer port is disconnected from the canister, make sure water does not enter the connector.
- Always clamp the transfer port tube and the canister tube before disconnecting them.
- Do not take the pump apart.
- Do not modify pump, canister, tubes, film or foam as any modification may significantly compromise the ability of Avance Solo Adapt NPWT System to deliver therapy.
- The Avance Solo Adapt Pump is not intended for use aboard aircraft. During air travel pause the pump and remove the batteries. Make sure that the dressing is not left without applied suction longer than the time determined by your healthcare professional.
- The potential for electromagnetic interference in all environments cannot be eliminated. Use caution if the pump is near electronic equipment such as anti-theft equipment or metal detectors and ensure proper function according to section 7 Handling instructions.

5. Description of Avance Solo Adapt NPWT System

Avance Solo Adapt NPWT System consists of Avance Solo Adapt Pump, Avance Solo Canister 50 ml, Avance Solo Adapt Foam, Avance Solo Adapt Film and Avance Solo Adapt Transfer Port.

The foam is placed in the wound, the film is then applied, the transfer port is applied on the film and then connected to the canister tube via the green connectors. The canister is attached to the pump. When the pump is started a negative pressure will be created, providing suction to the wound. Fluid from the wound will be transported to and collected in the canister. If the canister becomes full, you or your health care specialist will be able to change the canister, see Section 7.3 for instructions.

The pump is operated with a single push button and is battery powered. If the battery power becomes low, you or your health care specialist will be able to change the batteries, see Section 7.4 for instructions. The pump has auditory (beeps) and visual (indicator lights) notifications and alarms to let you know when therapy is being applied as intended or if there is an issue: Make sure that you place the pump so that you can detect the notifications and alarms. See Section 7 and 12 for more instructions and troubleshooting guides.

6. Daily life while on therapy with Avance Solo Adapt NPWT System

Can you move around while on the therapy? Based on your health condition, you should be able to move around and maintain your daily activities. Follow the instructions provided by your healthcare professional.

Will it be painful? When the dressing is first applied, and the pump is started, you may experience a slight pulling or drawing sensation from dressing contraction. If you experience any pain during therapy, please consult your healthcare professional for advice.

How often will the dressing be changed? How often your dressing needs to be changed depends on type of wound and how much fluid is collected from the wound. A healthcare professional will determine and give you information on how often your dressing will be changed. The dressing must be changed by a healthcare professional. A typical change frequency is 48 to 72 hours but no less than 3 times a week, or as instructed by healthcare professional.


When you are resting or asleep, place the pump in a safe position where it cannot be pulled off a table or cabinet onto the floor. Make sure that all tubes are placed in a position that minimizes risk of entrapment or strangulation.

Can I shower? Light showering is allowed but keep the pump away from water. If the pump is accidentally wet, disconnect the pump and canister and contact your healthcare professional. The dressing is water resistant but should not be exposed to jets of water.

For light showering, pause therapy by pressing and holding the green push button on the pump, and release after two (2) seconds. Clamp the canister tube and transfer port tube by positioning the slide clamps next to the green connectors and slide across the tubes until secured. Disconnect the canister tube from the transfer port tube. Ensure that the transfer port tube is not placed into contact with water.

Can I clean the pump? You can clean the pump by wiping with a damp cloth or with a non-abrasive detergent. Do not put the pump under running water.


7. Handling instructions

7.1. How do I know that Avance Solo Adapt NPWT System gives therapy as intended?

When Avance Solo Adapt NPWT System gives therapy as intended, the green push button on the pump will flash 2 times every minute and the dressing will have a wrinkled appearance and be firm to the touch. Regularly check that negative pressure is active by monitoring visual and audible notifications and alarms from the pump.

NOTE: When the pump is first started, the green push button on the pump flashes once every second for 15 minutes.

7.2. How do I know that the pump is paused?

When the pump is paused you will hear an audible notification as two short beeps and the green push button will stop flashing, all indicator lights will be off. The pump will repeat the two short beeps every 15 minutes as long as the pump is paused.

NOTE: If you do not restart the pump manually, it will automatically restart after 60 minutes.


7.3. How do I know when to change the canister?

The canister should be changed when you can see that it is full through the transparent window on the back of the canister. A Blockage alarm from the pump may also indicate that you need to change the canister. If no change is required, the canister may be used up to 14 days.


The **BLOCKAGE** indicator light flashes once every second, the pump repeatedly alarms with a beep and will then pause.

NOTE: A Blockage alarm may also be triggered by blockage in the tubes. Always make sure that the tubes are not clamped or kinked.


To replace the canister, perform the following steps


- If the pump is active, pause the pump by pressing down the green push button, and release after (2) seconds.
- Clamp the canister tube and transfer port tube by positioning the slide clamps next to the green connectors and slide across the tubes until secured. Clamping the tubes minimizes fluid leakage when you disconnect the canister tube from the transfer port tube.
- Disconnect the canister tube from the transfer port tube by squeezing the connector from both sides and pull apart.

- Remove the canister by pushing the spring buttons on both sides and pull.
- Attach a new canister to the pump by pushing the canister until it clicks on both sides and is locked in place.
- Be sure to attach the canister tube onto the back of the pump.

- To continue therapy, connect the canister tube to the transfer port tube.
- Make sure that the clamps on the canister and transfer port tubes are released. Restart the pump by pressing down the green push button, and release after two (2) seconds.
- Monitor that the negative pressure is active, the dressing should be contracted and firm to the touch.

7.4. How do I know when to change batteries?


The Avance Solo Adapt Pump is battery powered and will require you or responsible healthcare professional to change batteries, normally after 7 days or when the pump alarms for Low battery. When the battery power is reaching its end, the pump will indicate the following:

When 24 hours of battery power remain, the indicator light for **LOW BATTERY** flashes once every five (5) seconds.


When less than 4 hours of battery power remain, the indicator light for **LOW BATTERY** flashes once every second and the pump repeatedly alarms with a beep.

To replace batteries, perform the following steps

Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 13.

- If the pump is still active, pause the pump by pressing down the green push button, and release after two (2) seconds.
- Open the battery compartment on the back of the pump by sliding the lid. Remove the batteries. Insert new batteries, ensuring that the positive terminal (marked +) and negative terminal (marked -) of each battery matches the +/- label in the battery compartment. Close the battery compartment lid.
- When batteries are correctly inserted, the pump will notify with three beeps at different frequency tones: one beep with high, one beep with medium followed by one beep with low frequency tone. This is an automatic self-check that confirms that the batteries are correctly inserted in the pump and that the pump is ready to use.
- The pump will then enter pause mode, all indicator lights will be off. Restart the pump by pressing down the green push button, and release after two (2) seconds. Make sure that the negative pressure is active, the dressing should be contracted and firm to the touch.

7.5. How do I know that the therapy time of 14 days is complete?

The Avance Solo Adapt Pump is for single patient use, is battery powered and has a 14-day lifespan. When the therapy time of 14 days is reached, the pump will indicate the following:

All indicator lights and the green push button on the pump will flash with high intensity and the pump will notify with three beeps at different frequency tones: one beep with high, one beep with medium followed by one beep with low frequency tone.

7.6. What do I do if the green push button is accidentally pressed?

When the green push button is accidentally pressed, the will pump will notify with a beep. No action is required.

8. Disposal

When you have changed the batteries, dispose the batteries so that they can be recycled in accordance with requirements by local regulations, relevant state laws and the Waste Electrical and Electronic Equipment Directive (WEEE). When you have changed the canister, make sure that you dispose it as clinical waste in accordance with local regulations. Ask your healthcare professional for more information if you are unsure as to the safe disposal. You can also find more information about safe disposal at www.molnlycke.com/wastehandling.

9. Caution

Avance Solo Adapt NPWT System must be used in accordance with instructions provided in this Patient and Caregiver User Manual. Read these instructions before using the system and have them available during use. Failure to read and understand these instructions may lead to misuse of the system and improper performance. These instructions are a general guide for the use of the product. Specific medical situations must be addressed by a healthcare professional.

10. Material content

Film: polyethylene, polyurethane, soft silicone, polyacrylate adhesive

Foam: polyurethane

Canister: polycarbonate, polyurethane

Pump: polycarbonate, acrylonitrile butadiene styrene, thermoplastic elastomer

Transfer Port and tubes with clamp: Polyurethane, polyester, polyacrylate adhesive, polyolefin based thermoplastic elastomer, polyethylene

Connectors: acrylonitrile butadiene styrene copolymer, thermoplastic olefin, polyethylene

11. Other information

If any serious incident has occurred in relation to the use of Avance Solo Adapt NPWT System, it should be reported to Mölnlycke Health Care and to your local competent authority.

12. Troubleshooting

During therapy with Avance Solo Adapt NPWT System, it is important that you are aware of the audible and visual alarms and notifications displayed by the pump. This section explains the audible and visual alarms and notifications and gives guidance on how to troubleshoot and when to contact a healthcare professional.

AUDIBLE AND VISUAL DISPLAY		POSSIBLE CAUSE	HOW TO TROUBLESHOOT
	The LEAKAGE indicator light flashes once every second (therapy is still active). If there is extensive leak, the pump repeatedly alarms with a beep and will then pause (therapy is not active).		
	The BLOCKAGE indicator light flashes once every second (therapy is still active). The pump repeatedly alarms with a beep and will then pause.	The canister is full or there is a blockage in the tubes.	To correct a blockage: Ensure that the tubes are not clamped or kinked. If the canister is full, change the canister according to instructions in section 7.3. Press and hold the green start button for two (2) seconds to restart the pump. If blockage is not resolved, the pump will gain alarm for blockage and then pause therapy. If this happens, contact your healthcare professional.
	The LOW BATTERY indicator light flashes once every five (5) seconds. When less than 4 hours of battery power remain the LOW BATTERY indicator light flashes once every second and the pump repeatedly alarms with a beep.	The batteries are low.	To change batteries: Change the batteries according to instructions in section 7.4. Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 13. Press and hold the green start button for two (2) seconds to restart the pump.
	All indicator lights flash at the same time, once every second, and the pump repeatedly alarms with a beep.	Alarm for INTERNAL FAULT and the pump cannot be started.	Contact your healthcare professional or Mölnlycke Health Care.

13. Specifications Avance Solo Adapt Pump

Parameter	Specification
Nominal negative pressure	-125 mmHg
Maximum negative pressure	-150 mmHg
Mode of Operation	Continuous
Dimensions	Avance Solo Adapt Pump and Avance Solo Canister 50 ml 125x68x30 mm
Weight	Avance Solo Adapt Pump and Avance Solo Canister 50 ml < 130 g
Applied Part	Dressing, type BF
Battery	2 x AA 1.5V Energizer L91
IP22	Ingress protection effective against fingers and similar objects. Protected against dripping water when tilted at 15°. Classification only valid when battery lid is closed.
Storage	Temperature 5°C/41 F to 25°C/77F, ambient humidity 10% to 75% non-condensing, ambient pressure 700 hPa to 1060 hPa
Transport	Temperature -35°C/-31 F to 63°C/145 F, ambient humidity 10% to 90% non-condensing, ambient pressure 700 hPa to 1060 hPa
Operation	Temperature 5°C/41 F to 40°C/104 F, ambient humidity 15% to 90%, non-condensing, ambient pressure 700 hPa to 1060 hPa
Low priority alarm signal, Alarm Volume 60 dBA	Leakage Alarm, Blockage Alarm, Low Battery Alarm, Internal Fault Alarm.
Information signals with lower priority than alarm signals	Pause mode, Therapy mode, Invalid button press, Pump Self-check, End of therapy, Leakage, Blockage, Low Battery.
Essential Performance	Activation of low priority alarms within two hours if degradation of Nominal negative pressure. Negative pressure not exceeding Maximum negative pressure longer than five minutes.

14. Safety

Avance Solo Adapt NPWT System complies with the General Requirements for Safety of Medical Electrical Equipment (IEC 60601-1). Avance Solo Adapt NPWT System is intended for home care use (IEC 60601-1-11).

15. Electromagnetic compatibility

The Avance Solo Adapt Pump has been tested in accordance with the requirements of IEC 60601-1-2. Exceeding test levels can cause degradation of negative pressure or negative pressure exceeding specifications. The pump may fail to provide alarm signals.

The Avance Solo Adapt Pump is tested for use in Professional healthcare facility environment and Home healthcare environment.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) from the Avance Solo Adapt Pump. Otherwise, degradation of the performance of this equipment could result.

Avance® is a registered trademark of Mölnlycke Health Care AB

Toll free number:

USA 1-800-882-4582

Canada 1-800-494-5134

Patent: www.molnlycke.com/virtual-patent-marking/

UK Responsible Person:

Mölnlycke Health Care Ltd.,

Medlock Street, Oldham, Lancashire, OL1 3HS, UK

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