HEALTHCARE PROFESSIONAL INSTRUCTIONS FOR USE

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

www.molnlycke.com

www.molnlycke.com/symbols

STERILE EO	The device is sterilized using ethylene oxide
	Do not reuse
	Caution, see instructions for use
	Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 8.
	Follow instructions for use
	Temperature limitation
	Keep dry Keep away from rain
	Keep away from sunlight Keep away from heat
	Catalogue number
	Use by date / Expiry date
	Batch code
	Serial number

7.2. Avance Solo NPWT System – Normal use

Avance Solo Pump displays the following visual and audible signals to inform the user that the Avance Solo NPWT system is performing under normal operation.

AUDIBLE AND VISUAL DISPLAY	DESCRIPTION	OPERATION	COMMENT
	1. The pump activates for a short period. 2. All light indicators on the pump flash in turn. 3. The pump emits audible notifications: medium frequency followed by high frequency	Automatic self-check	Automatic self-check is performed when the batteries are correctly inserted in the pump and confirms that the pump is ready for use.
	1. At start-up of the pump, the green start button flashes once every second for 15 minutes. 2. During normal operation the green start button flashes two times every minute.	Therapy mode	1. To confirm that correct negative pressure is reached. 2. To confirm that the pump is operating correctly, and that therapy is maintained.
	1. The pump emits two short audible notifications 2. The audible notifications are repeated every 15 minutes as long as the therapy is paused.	Pause mode	To confirm that the pump and therapy have been paused.
	1. All light indicators flash with high intensity. 2. The pump emits three audible notifications: one high, one medium followed by one low frequency tone.	End of therapy	The therapy time of 14 days is complete.
	The pump emits short audible notification	Invalid button press	In case of an invalid button press

7.3. Avance Solo NPWT System – Alarms and trouble shooting

Avance Solo Pump displays the following visual signals and audible alarms to inform the user of risk for loss of therapy.

AUDIBLE AND VISUAL DISPLAY	DESCRIPTION	Possible Cause	Comments / Trouble Shooting
	1. The indication light for leakage flashes once per second. 2. If negative pressure is not maintained: • The light indicator for air leakage stays active, flashing once every second • The pump repeatedly emits an audible alarm	Leakage Alarm Negative pressure is not being established due to a leak.	To correct a leak, perform one or more of the following actions: Press around the dressing border to improve contact with the skin or if required add additional fixation strips around dressing borders. Ensure canister is securely attached to the pump. Ensure that the tube is securely attached to the canister. Ensure that the dressing tube is securely connected to the canister tube.
	1. The indication light for blockage flashes once every second. 2. The blockage condition persists: • The light indicator for blockage stays active, flashing once every second • The pump repeatedly emits an audible alarm	Blockage Alarm Negative pressure is not being established due to blockage	To correct a blockage, perform one or more of the following actions: Ensure that the tubes are not clamped. Ensure that the tubes are not full. If the canister is full perform a canister change according to instructions for canister change.
	1. The indication light for battery flashes once every five (5) seconds when up to 24 hours of battery time remain. 2. When less than 4 hours of battery time remain: • The light indicator for battery stays active, flashing once every second • The pump repeatedly emits an audible alarm	Low Battery Alarm	To correct a low battery alarm: Change the batteries according to instructions for battery change. Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 8.
	1. The indication lights for leakage, battery and blockage flash simultaneously once every second. 2. The pump repeatedly emits an audible alarm.	Internal Failure Alarm	The pump has an internal fault and cannot be started. Contact your healthcare professional or Mölnlycke Health Care.

8. Specifications Avance Solo Pump

Nominal negative pressure	-125 mmHg
Maximum negative pressure	-150 mmHg
Mode of operation	Continuous
Dimensions	Avance Solo Pump and Canister 50 ml 125x68x30 mm
Weight	Avance Solo Pump and Canister 50 ml < 130 g
Applied Part	Dressing, type BF
Battery	2x AA 1.5V Energizer L91
IP22	Ingress protection effective against fingers and similar objects. Protected against dripping water when tilted at 15°. Classification only valid when battery lid is closed.
Storage	Temperature 5°C/41 F to 25°C/77F, ambient humidity 10% to 75% non-condensing, ambient pressure 700 hPa to 1060 hPa
Transport	Temperature -35°C/-31 F to 63°C/145 F, ambient humidity 10% to 90% non-condensing, ambient pressure 700 hPa to 1060 hPa
Operation	Temperature 5°C/41 F to 40°C/104 F, ambient humidity 15% to 90%, non-condensing, ambient pressure 700 hPa to 1060 hPa
Low priority alarm signal, Alarm Volume 60 dB	Leakage Alarm, Blockage Alarm, Low Battery Alarm, Internal Failure Alarm.
Information signals with lower priority than alarm signals	Pause mode, Therapy mode, Invalid button press, Pump Self-check, End of therapy, Leakage, Blockage, Low Battery.
Essential Performance	Activation of low priority alarms within two hours if degradation of Nominal negative pressure. Negative pressure not exceeding Maximum negative pressure longer than five minutes.

9. Safety

Avance Solo NPWT System complies with the General Requirements for Safety of Medical Electrical Equipment (IEC 60601-1). Avance Solo NPWT System is intended for home care use (IEC 60601-1-11).

Conforms to CSA Std. ES60601-1, AAMI IEC Std. 60601-1-8 and AAMI Std. HA60601-1-11. Certified to CSA Std. C22.2 No. 60601-1-6, CSA Std. C22.2 No. 60601-1-8 and CSA Std. C22.2 No. 60601-1-11.

10. Electromagnetic compatibility

Avance Solo Pump is tested in accordance with the requirements of IEC 60601-1-2. Exceeding test levels can cause degradation of negative pressure or negative pressure exceeding specifications. The pump may fail to provide alarm signals.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to Avance Solo Pump. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic Emissions
Avance Solo Pump is intended for use in the electromagnetic environment specified below.

EMISSION TESTS	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	Avance Solo Pump uses RF energy only for its internal function
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	The RF emissions characteristic of the Avance Solo Pump make it suitable for use in home use environments and in hospitals, except for near high frequency surgical equipment

Electromagnetic Immunity

Avance Solo Pump is intended for use in the electromagnetic environment specified below.

IMMUNITY TESTS	BASIC EMC STANDARD OR TEST METHOD	IMMUNITY TEST LEVELS
IMMUNITY TESTS		Professional healthcare facility environment Home healthcare environment
Electrostatic Discharge	61000-4-2	± 8 kV contact ± 2 kV, ± 8 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	61000-4-3	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz
Proximity fields from RF wireless communications equipment	61000-4-3	30 cm minimum separation distance from radio transmitter
Rated power frequency magnetic fields	61000-4-8	30 A/m 50 Hz or 60 Hz

Recommended separation distances between RF transmitters and Avance Solo Pump IEC 60601-1-2:2014 Table 9

Frequency band (MHz)	Service	Maximum Power (W)	Minimum distance	Immunity Test Level (V/m)
380-390	TETRA 400	1.8	0.3 meter	27
430-470	GMRS 460, FRS 460	2	0.3 meter	28
704-787	TDE Band 13, 17	0.2	0.3 meter	9
800-960	GSM, TETRA, iDEN, CDMA, LTE	2	0.3 meter	28
1700-1990	GSM, CDMA, DECT, LTE, UMTS	2	0.3 meter	28
2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID, LTE	2	0.3 meter	28
5100-5800	WLAN 802.11 a/n	0.2	0.3 meter	9

11. Caution

Avance Solo NPWT System must be used in accordance with these Instructions for Use. Read these instructions before using the system and have them available during use. Failure to read and understand these instructions may lead to misuse of the system and/or poor performance. These instructions are a general guide for the use of the product. Specific medical situations must be addressed by a clinician.

12. Other information

If any serious incident has occurred in relation to Avance Solo NPWT System, it should be reported to Mölnlycke Health Care and to your local competent authority.

13. UK Responsible Person:

Australian sponsor address : Mölnlycke Health Care Pty. Ltd

USA : 1-800-882-4582

Medlock Street, Oldham, Lancashire, OL1 3HS, UK

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